



A station which enables you to make a fast and objective quality control on your samples...

# OPTO - AgriMetric

## WHEAT ANALYSIS STATION by image processing


Percentage of normal grains and shrivelled grains

Percentage of «yellow berry» on normal grains  
(non fully vitreous kernels)

Kernel biometrics

Broken kernel rate

Measurement of Thousand-Seed Weight



Identification of normal grains, shrivelled grains and broken grains

Grain to grain objective biometrics

Results available in a few seconds



**Industrial  
and  
ergonomic  
system!**

**Patented  
technology**

Wiser for its several years experience in seeds and kernels image analysis (rice, wheat, corn...), Optomachines team proposes various systems and units for defects detection, biometrics and TSW measurement.

For more details,  
please contact us.

# OPTOmachines

CONSEIL - CONCEPTION - FABRICATION  
DE SYSTÈMES DE MESURE ET DE CONTRÔLE

2 rue Richard Wagner

Tél./Fax : +(33) 4 - 73 - 64 - 10 - 06

63200 Riom - France

Mail: [optomachines@optomachines.fr](mailto:optomachines@optomachines.fr)

<http://www.optomachines.fr>

Dimensions (h x w x d):  
91 x 44 x 42 cm

Included:

- Class 2 balance.
- PC with 19 inches TFT screen.

